

Pasny



microneedles




Ultra-fine With Thinner Walls

Less Puncture Resistance

Painstaking manufacturing with the utmost attention to detail based on the very latest technology has resulted in this ultra-fine needle with thinner walls, which offer numerous benefits to both doctor and patient.

- Reduces patient pain and discomfort to the absolute minimum
- Reduces bleeding to the minimum
- Minimal inflammation, minimising patient recovery time
- Facilitates penetration of active ingredients





BENEFITS OF PASNY MICRONEEDLES

Very Sharp Tips

Ultra-thin needles that minimise discomfort and pain during insertion. Manufactured with painstaking sharpening of the tip and a special polishing system, they offer virtually no friction resistance.

Ultra-Thin Walls

Pasny and Paskin needles have ultra-thin walls, facilitating the injection of highly viscous products, allowing the use of smaller gauge needles.

Less Discomfort

The patient's sensation of pain is reduced by up to 58% compared to normal needles.

Reference	Needle	Needle length	Presentation
PS-016-01	34G	4.0 mm	100 u.
PS-016-01		8.0 mm	100 u.

Comparative with other marks

The thinner of the market

JBP nanoneedle34G	PASNY 34 G	
0,20	0,18	
0,10	0,10	

RESULTS AFTER A FACIAL MESOTHERAPY SESSION



Immediately after



Immediately after

